

## Practical procedures:

### ***Making a standard solution:***

**Stage 1:** transfers known mass of solid

- a) Weigh the sample bottle containing the solid on a (2 dp) balance
- b) Transfer to beaker\* and reweigh sample bottle
- c) Record the difference in mass

**Stage 2:** Dissolves in water

- a) Add distilled / deionised water
- b) Stir (with a glass rod) or swirl
- c) Until all solid has dissolved

**Stage 3:** Transfer, washing and agitation

- a) Transfer to volumetric / graduated flask. Allow if a clear description/diagram given eg long necked flask with 250 cm<sup>3</sup> mark
- b) With washings
- c) Make up to 250 cm<sup>3</sup> / mark with water
- d) Shakes/inverts/mixes

### ***Recrystallisation:***

Dissolve the product in the **minimum** volume of Hot water / solvent

Allow the solution to cool and allow crystals to form.

Filter off the pure product under reduced pressure / using a Buchner funnel and side arm flask

The crystals were compressed in the funnel

A little cold water was poured through the crystals

### ***Melting point determination***

Measure the melting point using melting point apparatus

Sharp melting point / melting point matches data source value

### ***Experiment to determine the solubility:***

Take a known volume of the saturated solution

Evaporate the filtrate to dryness

Weigh the residue

### ***Measuring temperature drop (eg in an enthalpy calculation):***

Place thermometer in an insulated cup.

Start a clock

Record the temperature every subsequent minute for 3 minutes

Add the **salt / solution** at minute 4

Record the temperature every subsequent minute from minute 5

Plot a graph of temperature vs time

Extrapolate back to time of mixing = 4 and determine the temperature

### **Measuring temperature rise (eg in an enthalpy calculation):**

Start a clock when **salt** is added to water

Start a clock when **salt** is added to water

Record the temperature every subsequent minute for about 5 minutes

Plot a graph of temperature vs time

Extrapolate back to time of mixing = 0 and determine the temperature

### **Preparing a pure sample of AgBr from a mixture of NaCl and NaBr**

#### **Stage 1: formation of precipitates**

- Add silver nitrate
- to form precipitates of AgCl and AgBr
- $\text{AgNO}_3 + \text{NaCl} \rightarrow \text{AgCl} + \text{NaNO}_3$
- $\text{AgNO}_3 + \text{NaBr} \rightarrow \text{AgBr} + \text{NaNO}_3$

#### **Stage 2: selective dissolving of AgCl**

- Add excess of dilute ammonia to the mixture of precipitates
- the silver chloride precipitate dissolves
- $\text{AgCl} + 2\text{NH}_3 \rightarrow \text{Ag}(\text{NH}_3)_2^+ + \text{Cl}^-$

#### **Stage 3: separation and purification of AgBr**

- Filter off the remaining silver bromide precipitate
- Wash to remove soluble compounds
- Dry to remove water

### **Producing a pH curve**

- Measure pH of 25cm<sup>3</sup> of the acid
- Add alkali in 1cm<sup>3</sup> portions
- Stir mixture
- Measure pH (after each addition)
- Repeat until 24cm<sup>3</sup>
- Add 0.5cm<sup>3</sup> near endpoint 24 – 26cm<sup>3</sup>
- Add alkali in 1cm<sup>3</sup> portions until in excess

### **Thin layer chromatography:**

#### **Set up:**

- Wearing plastic gloves to hold a TLC plate, draw a pencil line 1.5 cm from the bottom of the plate.
- Use a capillary tube to apply a very small drop of the solution of amino acids to the mid-point of the pencil line.
- Allow the spot to dry completely.
- Add the developing solvent to a depth of not more than 1 cm.
- Place your TLC plate in the developing tank.
- Allow the developing solvent to rise up the plate to the top.
- Remove the plate and quickly mark the position of the solvent front with a pencil.
- Allow the plate to dry in a fume cupboard.

#### **Developing:**

- Spray with developing agent
- Measure distances from initial pencil line to the spots (x)
- Measure distance from initial pencil line to solvent front line (y)
- $R_f \text{ value} = x / y$
- Amino acids have different polarities
- Therefore, have different retention on the stationary phase or different solubility in the developing solvent

### **How would you determine a suitable solvent for recrystallisation:**

Place a small amount of pure crystals in a test tube.

Add small volume (0.5–1.0 cm<sup>3</sup>) of possible solvent.

Shake/stir.

If pure crystals dissolve then solvent unsuitable.

If the pure crystals are insoluble in cold solvent then heat in a water bath.

If pure crystals do not dissolve = unsuitable.

If pure crystals dissolve partially add more solvent.

If/when completely dissolved place in ice water bath.

If crystals form = suitable.

Effectiveness determined by measuring melting point.

Purity indicated by melting point being sharp and close to data book melting point

### **Colorimetry to measure rates:**

#### **Setting up**

Calibrate the colorimeter with known concentrations

Plot a graph of absorbance against concentration

#### **Carrying out**

Fill a cuvette with the solution at regular time intervals during the reaction

Read the absorbance from the colorimeter

Read from the graph to find the concentration of the solution

Plot a graph of concentration against time


If it is a straight line = 0 order

Constant half lives = 1<sup>st</sup> order


Increasing half lives = 2<sup>nd</sup> order

### **Practical techniques:**


#### **Fractional distillation**




#### **Distillation**




#### **Reflux**



#### **Filtration**



#### **Filtration under reduced pressure**



#### **Gas collection**

